## OBTAIN IMPRESSIVE IMAGES

Match RiX70 DC high frequency X-Ray generator with any of the Trident image receptors and get more than you could have ever imagined.



### WHAT'S THE SENSE OF DC?

DC means amazing results with all the advantages of the latest electronic technology: rather than using alternating current, RiX70 DC applies a nearly constant potential to the tube. This constant potential supplies the tube with DC current, resulting in a constant stream of consistent radiation that maintains the preset kV values even with mains fluctuations.




SYSTEM SUPPLY		
Supply Voltage	AC100-240V	
Frequency	50/60 Hz	
Electric current	2.7 A @ 240V , 6.4 A @ 100V	
Max Power Consumption	640VA (3 AT fuse)	

X-RAY GENERATOR		
Tube Voltage	60 kVρ ± 8% - 65 kVρ ± 8% - 70 kVρ ± 8%	
Frequency	75 kHz	
Anodic current	7 mA± 10%	
Output radiation linearity	< 0.2 according to IEC 60601-2-63	
Duty cycle	1:30	
Nominal power	490 W (70 kVρ - 7 mA)	
Total filtration	1.5 mm Al eq. @ 70 kVp	
Focal spot	0.4 mm (IEC 60336)	
Leakage radiation at 1 m	<0,1 mGy/h	
Max thermal capacity of monoblock	130J (convection oil bath cooling)	
Distance focus to skin	20 cm	
Beam diameter X Ray at 200 mm	60mm (optional 36 x 45 mm rectangular collimator)	

X-KAY TUBE			
Canon	Skan X - CEI (Italy)		
D-045	0X70-4		
< 1 mm Al eq.	0.5 mm Al eq.		
Max. 3 A	1.5 – 2.1 A		
0.4 IEC 60336:2005			
Tungsten			
70 kV			
	Canon D-045 < 1 mm Al eq. Max. 3 A 0.4 IEC 60		

#### WALL MOUNTED VERSION



#### Trident S.r.l. Via Artigiani, 4 25014 Castenedolo Phone +39 030 2732485 info@trident-dental.com

# RiX70 DC High frequency intraoral X-Ray generator





#### **INSIDE RIX70 DC**



#### **AUTOMATIC VOLTAGE TRANSFORMER**

RiX70 DC appropriately works in a range of 100 to 240 V. The electrical system automatically detects the power wattage, allowing a safe and proper functioning with every voltage.



#### 0.4 MM FOCAL SPOT

The base for an accurate diagnostic is a good image. RiX70 DC with focal spot 0,4 mm gives you incredibly clear images with a high contrast range.



#### **REDUCTION OF RADIATION DOSE**

The high frequency generator allows radiation stability eliminating the soft radiation that does not produce images but can be harmful for the patient.



#### **TUBE VOLTAGE SELECTOR**

RiX70 DC operates either 60, 65 or 70 kV allowing to reach excellent results to satisfy your needs: Voltage 60 kV for caries detection or periodontal applications; 70 kV for endodontic applications, including visualization of canals, detection of root fractures and diagnosis of periapical lesions due to pulpal inflammation.



#### **AUTOMATIC SETTING OF EXPOSURE TIME**

From 0.01 to 2.0 seconds, depending on the patient and tooth characteristics

#### FOR BETTER RESULTS, YOU CAN EASILY:

- Select the exposure parameters according to the patient morphology, tooth type and image receptor.
- Adapt the radiation dose to lower or higher values.
- Modify the exposure time for digital or film.

## OPTIONAL



Install the remote timer outside the X-ray room. Conveniently set the exposure parameters, type of tooth, patient and image receptor and, safely shoot the X-ray.

Only available for the wall mounted version.

All local command devices are not necessary (control with coiled cord and the small keypad in the tube head support)



## MODERN DESIGN FOR PRACTICAL OPERATION

- Compact and light tubehead.
- Steady and strong support arm.
- Hand switch with 3 mt coiled cord.
- Additional keypad installed near the tube head, for manual adjustment of the exposure time and tooth type.
- Easy-to-read main control panel, with quickly identifiable icons.
- Easy and fast installation system.
- Available in wall mounted and mobile version.
- Length of the arm 60 and 80 cm.



